

BLOOD TRANSFUSION ARTEFACT RECORD	
NUMBER	AA-018
ITEM / TITLE	STERILE BLOOD ADMINISTRATION SET FOR BLOOD OR BLOOD COMPONENTS - BAXTER
PHYSICAL DESCRIPTION	Sealed half clear plastic / half paper bag containing a sterile blood administration set for blood or blood components with luer-lock adaptor. Blue writing on white paper side of the package includes information for use together with explanatory diagrams.
OTHER INFORMATION	<p>The package information states: BAXTER C2071 STERILE, NONPYROGENIC LOT 93 E 10U11 STERILISED 05 93 <i>Do not use if protectors are loose or missing.</i> <i>Discard after single patient use.</i> <i>Replacement is recommended within 24 hours.</i> <i>20 drops distilled water = 1 ml (approx.)</i> <i>15 drops blood = 1 ml (approx.)</i> <i>Directions for use: Use aseptic technique</i> <i>Cautions: 1. Maximum needle size for addition of medication through injection site: 23 SWG. 2. During pressure infusion ensure that both chambers remain full and air free.</i> <i>Manufactured by Baxter Healthcare Ltd. Thetford, Norfolk, England.</i> The label also lists 5 points for Retrograde Priming and 4 points for Normal Priming.</p>

PHOTOGRAPHS OF THE ARTEFACT	
	
	