



| BLOOD TRANSFUSION ARTEFACT RECORD | |
|-----------------------------------|---|
| NUMBER | AA-035 |
| ITEM / TITLE | A BOTTLE OF DRIED HUMAN ALBUMIN FOR TRANSFUSION |
| PHYSICAL DESCRIPTION | A 550ml Medical Research Council glass bottle of 'Dried Human Albumin for Transfusion with a silver metal cap cover and white label. Prepared and issued by the Blood Products Laboratory, Lister Institute, Elstree, Herts. |
| OTHER INFORMATION | <p>The label states:</p> <p><i>This bottle contains 24 grams of human albumin with 0.33 grams of Sodium Caprylate. Reconstitute by adding the volume of sterile pyrogen free distilled water, saline or glucose solution necessary to give the required protein concentration. Do not use vigorous shaking to hasten the dissolving of the solids.</i></p> <p><i>Reconstituted albumin MUST BE DISCARDED if not used within 3 hours.</i></p> <p>STORE IN A COOL DARK DRY PLACE</p> <p>Batch: EPA 619A/26</p> <p>Date: 27 SEP 1965</p> <p><i>Enter these details in patient's notes and quote in any enquiry.</i></p> <p><i>DO NOT USE if after adding fluid, a gel forms or solution is incomplete.</i></p> <p><i>Expiry Date: 27 SEP 1970</i></p> <p>Note: The 'DO NOT USE' section of the label is underlined.</p> |

| PHOTOGRAPHS OF THE ARTEFACT | |
|---|--|
|  |  |