

BLOOD TRANSFUSION ARTEFACT RECORD	
NUMBER	AA-030
ITEM / TITLE	A BOTTLE OF DRIED POOLED HUMAN PLASMA FOR TRANSFUSION
PHYSICAL DESCRIPTION	A 550ml Medical Research Council glass bottle of Dried Human Plasma with silver metal cap cover and white label. Stated to have been 'Prepared and Issued on behalf of the Department of Health & Social Security by the Blood Products Laboratory, Lister Institute, Elstree, Herts.'
OTHER INFORMATION	<p>The label states:</p> <p><i>Contains dried solids from 400 ml. unfiltered pooled citrate plasma. May be given to patients of ANY BLOOD GROUP. To reconstitute add 400 ml. pyrogen free distilled water.</i></p> <p><i>Protein 5% Citrate 1% Glucose 0.6%</i></p> <p><i>To administer with a piercing needle transfusion set, replace metal cap and rubber wad and SCREW CAP DOWN FIRMLY.</i></p> <p><i>Reconstituted plasma MUST BE DISCARDED if not used within 3 hours. STORE IN A COOL DRY DARK PLACE</i></p> <p><i>Issued by Regional Transfusion Centre.....</i></p> <p><i>Batch: 6D 4722/2</i></p> <p><i>Date: MAR 1970</i></p> <p><i>Enter these details in patient's notes and quote in any enquiry.</i></p> <p><i><u>DO NOT USE if after adding fluid, a gel forms or solution is delayed beyond 10 mins.</u></i></p> <p><i>Expiry: FEB 1978</i></p>

